



Specifications

	A8 Mini	A8 Pro
Model		
Dimensions	714mm(L) × 508mm(W) × 609mm(H)	750/1200mm(L) × 730mm(W) × 890mm(H) (Configurable chassis)
Weight	41.2kg	128kg
Power Supply	~100-240 VAC, 50/60Hz	
Arm Configurations	2-channel flexible pipetting channel module	
	8-channel fixed pipetting channel module	
	Plate transportation module	
Deck Capacity	11	20-35
Pipetting Principle	Air displacement	
Pipetting Range	1μL-1000μL	
Liquid Level Detection	Pressure liquid level detection	
Precision	2μL CV<5% ; 200μL CV<0.75%	
Precision of Arm Movement	± 0.1 mm on X,Y,Z axes	
Environmental Control	None. Can be placed in the biosafety cabinet	Optional UV & HEPA

* This data is estimated for reference only



A8 Series

Liquid Handling Workstation
Affordable Benchtop Solution

✦ Flexibility

✦ Safety

✦ Compact

✦ Traceability



FOLLOW US

[linkedin.com/company/nayo-biotec](https://www.linkedin.com/company/nayo-biotec)

[X.com/NayoBiotec](https://x.com/NayoBiotec)

facebook.com/nayobiotec

www.nayolab.com



By Precision, For Simplicity

www.nayolab.com

A8 Series


Liquid Handling Workstation

Modular and scalable design, adaptable to a wide range of flexible needs

- 8-channel and 2-channel pipetting modules, along with a Plate transportation module, can be flexibly configured
- Customizable workbench length to meet varying throughput demands
- Integration of various in-house and third-party modules to accommodate diverse application needs






8-channel fixed pipetting channel module



2-channel flexible pipetting channel module




Plate transportation module



Comprehensive process monitoring significantly extends walk-away time

- Feature pressure liquid level detection (pLLD) for real-time pipetting status monitoring
- Status monitoring is available during the transportation process with the plate transportation module



Graphical Interface, Easy to Use

- Simple process flow with one-click startup
- Intuitive drag-and-drop graphical interface for ease of learning



Integrated Functional Device



Plate Carousel




Screw Cap Capper/Decaper

Integrated Modules



Temperature Control Module



Magnetic Plate



Thermocycler Module

3rd Party Device






Plate Washer




Microplate Readers



Centrifuge



Refrigerator



Real Time Cycler

.....

Applications

Cell Biology

Molecular Biology

Toxicology

Drug Discovery

Clinical Diagnostics

Next Generation Sequencing

Synthetic Biology

Proteomics

Biobanking